

# DRAQ5™



Rapide soluzioni per la Ricerca Avanzata

rapid far-red emitting fluorescent DNA dye used for LIVE, permeabilized and fixed cell analysis

## DRAQ5 (5mM)

50µl	DR50050	250 Imaging Assays (slide); 100 Flow Cytometry Assays
200 µl	DR50200	1000 Imaging Assays (slide); 6600 HCS Assays (384-MTP); 400 Flow Cytometry Assays
1ml	DR51000	33000 HCS Assays (384-MTP)


## FEATURES AND BENEFITS

- Rapid staining of dsDNA/nuclei of LIVE or fixed cells
- Easy to use - no lyse, no wash, no RNase needed
- Ideal for use with GFP & FITC labels - DRAQ5 fluoresces in the far-red region
- Cost-effective analysis
- Fully compatible with CyGel™ thermoreversible gel BOS-CS20500

**Background information:** DRAQ5™ is a novel A-T intercalating anthraquinone. It has been chemically designed to enter the membranes of living cells, target dsDNA very specifically and contains a far-red fluorescing chromophore. DRAQ5™ is supplied in aqueous solution and is incredibly stable. It should be stored in stock solution at +4°C. **Major Applications.** DRAQ5™ has many applications and is highly compatible with existing protocols across a wide range of instrumentation platforms:

- nucleated cell gating in flow cytometry
- cell cycle and ploidy analysis / proliferation index
- nuclear counterstain and nuclear:cytoplasm segmentation in HCS

## DRAQ5™ SPECTRAL PROFILE



DRAQ5™ is a Trademark of Biostat Inc.

## Instrument Settings:

Suggested optimal emission filters include Cy5 (not recommended)

- for cell cycle analysis), 695LP, 730/50BP or 780/60BP.
- No compensation is needed with common FITC/PE antibody pairs in flow cytometry experiments.
- DRAQ5 allows multi-line imaging/cytometry when combined with UV & vis. range fluorochromes:
  - up to 4 colours in multi-colour microscopy
  - up to 6 colours using a 2-colour flow cytometer
- Photobleaching is extremely low allowing detection of weakly expressed reporter signals without loss of counterstain detail.
- Due to the excellent signal:noise there is no need to wash out DRAQ5 prior to imaging reducing pipetting steps and unnecessary sample processing.

## DRAQ5™ EXAMPLE PROTOCOLS

### Imaging

#### Cell staining with DRAQ5 for DNA cell cycle analysis

Last Updated: 03/09/2012

##### CELL STAINING WITH DRAQ5™ FOR DNA CELL CYCLE ANALYSIS BY FLOW CYTOMETRY OR BY CELL IMAGING -

Reagents required:

- DRAQ5™
- Phosphate Buffered Saline (PBS, without sodium azide) or other culture medium

1. Read the supplied Material Safety Data Sheet before handling DRAQ5™
2. Since no washing step is required, DRAQ5 will usually be the final staining procedure, after any cell treatment or labelling, prior to analysis.
3. Prepare cells for staining with DRAQ5™. Resuspend cells in appropriate buffer such as PBS at a concentration of  $\leq 4 \times 10^5$  / ml in a test tube. For adherent cells estimate the number of cells based on confluence level or tissue section dimensions.
4. Add DRAQ5™ directly as supplied following the 10  $\mu$ M or 20  $\mu$ M pipetting volumes in table 1. (For simple flow cytometric gating of nucleated cells the concentration of DRAQ5™ may be reduced). This will be as an overlay for adherent cells / tissue sections, added to the chamber/well liquid directly or in fresh medium following a wash step.
5. Gently mix and then incubate for 5-30 minutes at room temperature. nb: Protect from the light during this incubation period if other (immuno-) fluorescent stains have been applied to the cells, prior to the DRAQ5™ labelling, and which may otherwise suffer photo-bleaching. DRAQ5™ staining is accelerated at 37°C and maybe reduced to 1-3 min. DRAQ5™ stains intact, live, fixed, permeabilized and dead cells.
6. Cells can be analysed directly without further treatment or washing.

It is also important to consider the combinations of fluorochromes and filters for the experiment:

**EXCITATION:** DRAQ5™ may be excited by wavelengths from 488 nm (in flow cytometry) and up to 647 nm ( $Ex\lambda_{max}$  646 nm). Despite low absorbance at 488 nm this excitation may offer optimal CVs for flow cytometric cell cycle analysis whilst allowing convenient combination with FITC and R-PE conjugates and EGFP.

**EMISSION:** this starts at 665 nm ( $Em\lambda_{max}$  681 nm / 697 nm intercalated to dsDNA). Suitable filters include 695LP, 715LP or 780 LP. For cell cycle analysis it is recommended to choose a filter (such as 715 LP) which excludes a significant proportion of signal from the small fraction of unbound DRAQ5™.

Table 1: Ready reckoner for volumes of DRAQ5™ (5mM) required for various cell concentrations:

Cell sample preparation:		VOLUME OF DRAQ5™ (AS SUPPLIED) REQUIRED FOR A CONCENTRATION OF:		
No. of cells:	in volume:	5 $\mu$ M	10 $\mu$ M	20 $\mu$ M
$1 \times 10^6$	2500 $\mu$ l	2.5 $\mu$ l	5 $\mu$ l	10 $\mu$ l
$4 \times 10^5$	1000 $\mu$ l	1 $\mu$ l	2 $\mu$ l	4 $\mu$ l
$2 \times 10^5$	500 $\mu$ l	0.5 $\mu$ l	1 $\mu$ l	2 $\mu$ l
$1 \times 10^5$	250 $\mu$ l	0.25 $\mu$ l	0.5 $\mu$ l	1 $\mu$ l
$5 \times 10^4$	125 $\mu$ l	0.13 $\mu$ l	0.25 $\mu$ l	0.5 $\mu$ l

#### LIVE cell staining with DRAQ5™ for nuclear visualisation

Last Updated: 03/09/2012

##### LIVE CELL STAINING WITH DRAQ5™ FOR NUCLEAR VISUALIZATION BY HCS IMAGING PLATFORM OR LASER SCANNING CONFOCAL / EPIFLUORESCENCE MICROSCOPY

Reagents required:

- DRAQ5™
- Phosphate Buffered Saline (PBS, without sodium azide) or other culture medium

1. Read the supplied Material Safety Data Sheet before handling DRAQ5™
2. DRAQ5™ is usually added as the last stain in a labelling procedure since no washing is required or conveniently in assay medium for a live cell assay.
3. Prepare cells for staining with DRAQ5™. Resuspend cells in appropriate buffer such as PBS at a concentration of  $\leq 4 \times 10^5$  / ml in a test tube. For adherent cells estimate the number of cells based on confluence level or tissue section dimensions.
4. Add DRAQ5™ directly as supplied following the 5  $\mu$ M or 10  $\mu$ M pipetting volumes in table 1. This will be as an overlay for adherent cells / tissue sections, added to the chamber/well liquid directly or in fresh medium following a wash step.
5. Gently mix and then incubate for 5-30 minutes at room temperature. For time-lapsed assays (e.g. studying translocation of an EGFP tagged protein) DRAQ5™ may be added to the assay medium for the duration of the assay (typically 0.5 - 3 hr.) at 1  $\mu$ M prior to any agonist / antagonist addition.

nb: Protect from the light during this incubation period if other (immuno-) fluorescent stains have been applied to the cells, prior to the DRAQ5™ labelling, and which may otherwise suffer photo-bleaching. DRAQ5™ staining is accelerated at 37°C and maybe reduced to 1-3 min. DRAQ5™ stains intact, live, permeabilized and dead cells. Cells can be analysed directly without further treatment or washing. It is also important to consider the combinations of fluorochromes and filters for the experiment:

**EXCITATION:** DRAQ5™ may be sub-optimally excited by wavelengths from 488 nm (in flow cytometry) and up to 647 nm ( $Ex\lambda_{max}$  646 nm). Typically, for cell imaging, excitation is typically performed with either 561 nm, 633 nm or 647 nm wavelengths.

**EMISSION:**  $Em\lambda_{max}$  681 nm / 697 nm intercalated to dsDNA. Suitable filters include 695LP, 715LP or 780 LP. DRAQ5™ has no spectral emission overlap with FITC, R-PE, GFP and many other fluorescing proteins allowing image acquisition in one scan.

## Fixed cell/tissue staining with DRAQ5 for nuclear visualisation

Last Updated: 03/09/2012

### FIXED CELL/TISSUE STAINING WITH DRAQ5™ FOR NUCLEAR VISUALIZATION BY HCS IMAGING PLATFORM OR LASER SCANNING CONFOCAL / EPIFLUORESCENCE MICROSCOPY

Reagents required:

- DRAQ5™
- Phosphate Buffered Saline (PBS, without sodium azide) or other culture medium
- Paraformaldehyde

**A. SEPARATE FIXATION AND COUNTERSTAINING STEPS** (typically for cells / tissues where one or more external (immuno-) fluorescent stains will be applied):

1. Read the supplied Material Safety Data Sheet before handling DRAQ5™
2. DRAQ5 is usually added last in a procedure as no washing is required.
3. Prepare the sample for fixation and subsequent staining with DRAQ5.
4. Prepare separate working solutions of 4% formaldehyde and 5  $\mu$ M DRAQ5 in PBS, pipetting 1  $\mu$ l of DRAQ5, as supplied, into 1000  $\mu$ l of PBS.
5. Overlay the slide or chamber/well with the 4% formaldehyde solution. Incubate for 15-30 minutes at room temperature / 37°C.
6. Gently aspirate the formaldehyde solution, and wash with PBS.
7. Perform any necessary permeabilization, (immuno-)staining and blocking steps.
8. Overlay the washed, aspirated sample with 5 $\mu$ M DRAQ5 solution. Incubate for 15-30 minutes at room temperature / 37°C, in the dark.

NOTE: Protect from the light during these incubation periods if other (immuno-) fluorescent stains have been applied to the samples, which may otherwise suffer photo-bleaching. Samples can be analysed directly without further treatment or washing.

DRAQ5™ staining is accelerated at 37°C and incubation time may be reduced but this should be checked by titration and for each cell type. DRAQ5™ stains all nucleated cells.

It is also important to consider the combinations of fluorochromes and filters:

**EXCITATION:** DRAQ5™ may be sub-optimally excited by wavelengths from 488 nm (in flow cytometry) and up to 647 nm ( $Ex\lambda_{max}$  646 nm). Typically, for cell imaging, excitation is typically performed with either 561 nm, 633 nm or 647 nm wavelengths.

**EMISSION:**  $Em\lambda_{max}$  681 nm / 697 nm intercalated to dsDNA. Suitable filters include 695LP, 715LP or 780 LP. DRAQ5™ has no spectral emission overlap with FITC, R-PE, GFP and many other fluorescing proteins allowing image acquisition in one scan.

It should be possible to differentially segment the nucleus and cytoplasm.

**B. FOR A COMBINED FIXATION AND COUNTER-STAINING STEP** (typically for cells / tissues expressing an endogenous fluorescent protein as the only analyte e.g. translocation of a GFP-tagged transcription factor), simplifying the protocol:

1. Read the supplied Material Safety Data Sheet before handling DRAQ5™
2. DRAQ5™ is usually added last in a labeling as washing is not required.
3. Prepare the sample for fixation and subsequent staining with DRAQ5™.
4. Prepare separate working solutions of 8% formaldehyde (FA) and 10  $\mu$ M DRAQ5 in PBS, pipetting 2  $\mu$ l of DRAQ5, as supplied, into 1000  $\mu$ l of PBS.
5. Overlay the slide or chamber/well with equal volumes of formaldehyde and DRAQ5 solutions. Alternatively, prepare a pre-mix of the DRAQ5 and FA working solutions to simplify and speed the workflow for multiple samples. Overlay the cells with this mixture. Incubate for 15-30 minutes at room temperature / 37°C, in the dark.

NOTE: Protect from light during the incubation period if other (immuno-) fluorescent stains are present, which may otherwise suffer photo-bleaching. Samples can be analysed directly without further treatment or washing. DRAQ5™ staining is accelerated at 37°C and incubation time may be reduced but this should be checked by titration and for each cell type. DRAQ5™ stains all nucleated cells. It is important to consider

combinations of fluorochromes / filters for the experiment:

**EXCITATION:** DRAQ5™ may be sub-optimally excited by wavelengths from 488 nm (in flow cytometry) and up to 647 nm ( $Ex\lambda_{max}$  646 nm). Typically, for cell imaging, excitation is typically performed with either 561 nm, 633 nm or 647 nm wavelengths.

**EMISSION:**  $Em\lambda_{max}$  681 nm / 697 nm intercalated to dsDNA. Suitable filters include 695LP, 715LP or 780 LP. DRAQ5™ has no spectral emission overlap with FITC, R-PE, GFP and many other fluorescing proteins allowing image acquisition in one scan. It should be possible to differentially segment the nucleus and cytoplasm.

## Cytometry

### Cell staining with DRAQ5 for DNA cell cycle analysis

Last Updated: 03/09/2012

#### CELL STAINING WITH DRAQ5™ FOR DNA CELL CYCLE ANALYSIS BY FLOW CYTOMETRY OR BY CELL IMAGING

Reagents required:

- DRAQ5™
- Phosphate Buffered Saline (PBS, without sodium azide) or other culture medium

1. Read the supplied Material Safety Data Sheet before handling DRAQ5™
2. Since no washing step is required, DRAQ5 will usually be the final staining procedure, after any cell treatment or labelling, prior to analysis.
3. Prepare cells for staining with DRAQ5™. Resuspend cells in appropriate buffer such as PBS at a concentration of  $\leq 4 \times 10^5$  / ml in a test tube. For adherent cells estimate the number of cells based on confluence level or tissue section dimensions.
4. Add DRAQ5™ directly as supplied following the 10  $\mu$ M or 20  $\mu$ M pipetting volumes in table 1. (For simple flow cytometric gating of nucleated cells the concentration of DRAQ5™ may be reduced). This will be as an overlay for adherent cells / tissue sections, added to the chamber/well liquid directly or in fresh medium following a wash step.
5. Gently mix and then incubate for 5-30 minutes at room temperature.

nb: Protect from the light during this incubation period if other (immuno-) fluorescent stains have been applied to the cells, prior to the DRAQ5™ labelling, and which may otherwise suffer photo-bleaching.

DRAQ5™ staining is accelerated at 37°C and maybe reduced to 1-3 min. DRAQ5™ stains intact, live, fixed, permeabilized and dead cells. Cells can be analysed directly without further treatment or washing.

It is also important to consider the combinations of fluorochromes and filters for the experiment:

**EXCITATION:** DRAQ5™ may be excited by wavelengths from 488 nm (in flow cytometry) and up to 647 nm ( $Ex\lambda_{max}$  646 nm). Despite low absorbance at 488 nm this excitation may offer optimal CVs for flow cytometric cell cycle analysis whilst allowing convenient combination with FITC and R-PE conjugates and EGFP.

**EMISSION:** this starts at 665 nm ( $Em\lambda_{max}$  681 nm / 697 nm intercalated to dsDNA). Suitable filters include 695LP, 715LP or 780 LP. For cell cycle analysis it is recommended to choose a filter (such as 715 LP) which excludes a significant proportion of signal from the small fraction of unbound DRAQ5™.

#### DRAQ5™ is intended for research purposes only.

1. Detection of cytoplasmic nucleophosmin expression by imaging flow cytometry. Grimwade L, Gudgin E, Bloxham D, Bottley G, Vassiliou G, Huntly B, Scott MA, Erber WN. *Cytometry A*. 2012 Oct;81(10):896-900. doi: 10.1002/cyto.a.22116. Epub 2012 Sep 11. PMID: 22968966
2. Rescue from a two hit, high-throughput model of neurodegeneration with N-acetyl cysteine. Unnithan AS, Choi HJ, Tittler AM, Posimo JM, Leak RK. *Neurochem Int*. 2012 Aug;61(3):356-68. doi: 10.1016/j.neuint.2012.06.001. Epub 2012 Jun 9. PMID: 22691629
3. Mouse embryonic stem cell adherent cell differentiation and cytotoxicity assay. Barrier M, Chandler K, Jeffay S, Hoopes M, Knudsen T, Hunter S. *Methods Mol Biol*. 2012;889:181-95. doi: 10.1007/978-1-61779-867-2\_11. PMID: 22669665
4. A novel flow cytometric approach to distinguish circulating endothelial cells from endothelial microparticles: relevance for the evaluation of endothelial dysfunction. Lanuti P, Santilli F, Marchisio M, Pierdomenico L, Vitacolonna E, Santavenere E, Iacone A, Davi G, Romano M, Miscia S. *J Immunol Methods*. 2012 Jun 29;380(1-2):16-22. doi: 10.1016/j.jim.2012.03.007. Epub 2012 Mar 31. PMID: 22484509
5. Histone deacetylases 2 and 9 are coexpressed and nuclear localized in human molar odontoblasts in vivo. Klinz FJ, Korkmaz Y, Bloch W, Raab WH, Addicks K. *Histochem Cell Biol*. 2012 Feb 2. PMID: 22297573
6. Red/far-red fluorescing DNA-specific anthraquinones for nucl:cyto segmentation and viability reporting in cell-based assays. Edward R. *Methods Enzymol*. 2012;505:23-45. doi: 10.1016/B978-0-12-388448-0.00010-3. PMID: 22289446
7. Fluorescent microscopy of viable *Batrachochytrium dendrobatidis*. Friesen LR, Kuhn RE. *J Parasitol*. 2012 Jun;98(3):509-12. doi: 10.1645/GE-2973.1. Epub 2012 Jan 18. PMID: 22257116
8. Image Analysis of Colocalization of Nuclear DNA and GFP Labelled HIF-1 $\alpha$  in Stable Transformsants. Goto T, Sato M, Takahashi E, Nomura Y. *Curr Pharm Biotechnol*. 2012 Nov 1;13(14):2547-50. PMID: 22039803
9. Antifungal activity of *Coriandrum sativum* essential oil, its mode of action against *Candida* species and potential synergism with amphotericin B. Silva F, Ferreira S, Duarte A, Mendonça DI, Domingues FC. *Phytomedicine*. 2011 Dec 15;19(1):42-7. doi: 10.1016/j.phymed.2011.06.033. Epub 2011 Jul 23. PMID: 21788125
10. Evaluation of 309 environmental chemicals using a mouse embryonic stem cell adherent cell differentiation and cytotoxicity assay. Chandler KJ, Barrier M, Jeffay S, Nichols HP, Kleinstreuer NC, Singh AV, Reif DM, Sipes NS, Judson RS, Dix DJ, Kavlock R, Hunter ES 3rd, Knudsen TB. *PLoS One*. 2011;6(6):e18540. doi: 10.1371/journal.pone.0018540. Epub 2011 Jun 7. PMID: 21666745
11. Mouse embryonic stem cell adherent cell differentiation and cytotoxicity (ACDC) assay. Barrier M, Jeffay S, Nichols HP, Chandler KJ, Hoopes MR, Slentz-Kesler K, Hunter ES 3rd. *Reprod Toxicol*. 2011 May;31(4):383-91. doi: 10.1016/j.reprotox.2011.01.007. Epub 2011 Feb 3. PMID: 21296659
12. Short exposure to the DNA intercalator DRAQ5 dislocates the transcription machinery and induces cell death. Richard E, Causse S, Spriet C, Fourné N, Trinel D, Darzacq X, Vandebunder B, Heliot L. *Photochem Photobiol*. 2011 Jan-Feb;87(1):256-61. doi: 10.1111/j.1751-1097.2010.00852.x. Epub 2010 Dec 10. PMID: 21175643
13. Influence of the live cell DNA marker DRAQ5 on chromatin-associated processes. Mari PO, Verbiest V, Sabbioneda S, Gourdin AM, Wijgers N, Dinant C, Lehmann AR, Vermeulen W, Giglia-Mari G. *DNA Repair (Amst)*. 2010 Jul 1;9(7):848-55. doi: 10.1016/j.dnarep.2010.04.001. Epub 2010 May 2. PMID: 20439168
14. The use of DRAQ5 to monitor intracellular DNA in *Escherichia coli* by flow cytometry. Silva F, Lourenço O, Pina-Vaz C, Rodrigues AG, Queiroz JA, Domingues FC. *J Fluoresc*. 2010 Jul;20(4):907-14. doi: 10.1007/s10895-010-0636-y. Epub 2010 Mar 30. PMID: 20352307
15. Use of DNA-specific anthraquinone dyes to directly reveal cytoplasmic and nuclear boundaries in live and fixed cells. Edward R. *Mol Cells*. 2009 Apr 30;27(4):391-6. doi: 10.1007/s10059-009-0066-3. Epub 2009 Apr 13. Review. PMID: 19390818
16. Induction of DNA damage response by the supravital probes of nucleic acids. Zhao H, Traganos F, Dobrucki J, Wlodkowic D, Darzynkiewicz Z. *Cytometry A*. 2009 Jun;75(6):510-9. doi: 10.1002/cyto.a.20727. PMID: 19373929

17. Use of internally nuclease-protected single-strand DNA oligonucleotides and silencing of the mismatch repair protein, MSH2, enhances the replication of corrected cells following gene editing. Papaioannou I, Disteler P, Owen JS. *J Gene Med.* 2009 Mar;11(3):267-74. doi: 10.1002/jgm.1296. PMID: 19153972
18. Microglial cell migration stimulated by ATP and C5a involve distinct molecular mechanisms: quantification of migration by a novel near-infrared method. Miller AM, Stella N. *Glia.* 2009 Jun;57(8):875-83. doi: 10.1002/glia.20813. PMID: 19053059
19. Oxidative burst assessment and neutrophil-platelet complexes in unlysed whole blood. Avendaño A, Sales-Pardo I, Marin L, Marin P, Petriz J. *J Immunol Methods.* 2008 Dec 31;339(2):124-31. doi: 10.1016/j.jim.2008.09.003. Epub 2008 Sep 25. PMID: 18823984
20. Effect of proliferating cell nuclear antigen ubiquitination and chromatin structure on the dynamic properties of the Y-family DNA polymerases. Sabbioneda S, Gourdin AM, Green CM, Zotter A, Giglia-Mani G, Houtsmuller A, Vermeulen W, Lehmann AR. *Mol Biol Cell.* 2008 Dec;19(12):5193-202. doi: 10.1091/mbc.E08-07-0724. Epub 2008 Sep 17. PMID: 18799611
21. Hairy cell leukemia with unusual loss of CD103 in a subset of the neoplastic population: immunophenotypic and cell cycle analysis by flow cytometry. Yuan CM, Yang LJ. *Int J Clin Exp Pathol.* 2008 Jan 1;1(4):381-6. PMID: 18787608
22. DRAQ5 labeling of nuclear DNA in live and fixed cells. Smith PJ, Wiltshire M, Errington RJ. *Curr Protoc Cytom.* 2004 May;Chapter 7:Unit 7.25. doi: 10.1002/0471142956.cy0725s28. PMID: 18770802
23. Analysis of the piggyBac transposase reveals a functional nuclear targeting signal in the 94 c-terminal residues. Keith JH, Fraser TS, Fraser MJ Jr. *BMC Mol Biol.* 2008 Aug 11;9:72. doi: 10.1186/1471-2199-9-72. PMID: 18694511
24. Multispectral imaging of hematopoietic cells: where flow meets morphology. McGrath KE, Bushnell TP, Palis J. *J Immunol Methods.* 2008 Jul 31;336(2):91-7. doi: 10.1016/j.jim.2008.04.012. Epub 2008 May 15. Review. PMID: 18539294
25. Interaction of a DNA intercalator DRAQ5, and a minor groove binder SYTO17, with chromatin in live cells--influence on chromatin organization and histone-DNA interactions. Wojcik K, Dobrucki JW. *Cytometry A.* 2008 Jun;73(6):555-62. doi: 10.1002/cyto.a.20573. PMID: 18459157
26. DRAQ5-based, no-lyse, no-wash bone marrow aspirate evaluation by flow cytometry. Allan RW, Ansari-Lari MA, Jordan S. *Am J Clin Pathol.* 2008 May;129(5):706-13. doi: 10.1309/D2M6Q7BYHWAM1FYQ. PMID: 18426729
27. Total nucleated cell differential for blood and bone marrow using a single tube in a five-color flow cytometer. Björnsson S, Wahlström S, Norström E, Bernevi I, O'Neill U, Johansson E, Runström H, Simonsson P. *Cytometry B Clin Cytom.* 2008 Mar;74(2):91-103. PMID: 18061952
28. Baculovirus-mediated immediate-early gene expression and nuclear reorganization in human cells. Laakkonen JP, Kaikkonen MU, Ronkainen PH, Ihalainen TO, Niskanen EA, Hakkinen M, Salminen M, Kulomaa MS, Ylä-Herttua S, Airene KJ, Vihinen-Ranta M. *Cell Microbiol.* 2008 Mar;10(3):667-81. Epub 2007 Nov 27. PMID: 18042259
29. Simultaneous measurement of nucleated cell counts and cellular differentials in rat bone marrow examination using flow cytometer. Kurata M, Iidaka T, Hamada Y, Horii I. *J Toxicol Sci.* 2007 Aug;32(3):289-99. PMID: 17785944
30. Optimized isolation enables ex vivo analysis of microglia from various central nervous system regions. de Haas AH, Boddeke HW, Brouwer N, Biber K. *Glia.* 2007 Oct;55(13):1374-84. PMID: 17661344
31. DRAQ5: improved flow cytometric DNA content analysis and minimal residual disease detection in childhood malignancies. Swerts K, Van Roy N, Benoit Y, Laureys G, Philippe J. *Clin Chim Acta.* 2007 Apr;379(1-2):154-7. Epub 2006 Dec 20. PMID: 17229416
32. Use of flow cytometry to follow the physiological states of microorganisms in cider fermentation processes. Herrero M, Quirós C, García LA, Díaz M. *Appl Environ Microbiol.* 2006 Oct;72(10):6725-33. PMID: 17021224
33. Spectral analysis of the DNA targeting bisalkylaminoanthraquinone DRAQ5 in intact living cells. Njoh KL, Patterson LH, Zloh M, Wiltshire M, Fisher J, Chappell S, Ameer-Beg S, Bai Y, Matthews D, Errington RJ, Smith PJ. *Cytometry A.* 2006 Aug 1;69(8):805-14. PMID: 16969814
34. DNA labeling in living cells. Martin RM, Leonhardt H, Cardoso MC. *Cytometry A.* 2005 Sep;67(1):45-52. PMID: 16082711
35. Advanced microscopy solutions for monitoring the kinetics and dynamics of drug-DNA targeting in living cells. Errington RJ, Ameer-Beg SM, Vojnovic B, Patterson LH, Zloh M, Smith PJ. *Adv Drug Deliv Rev.* 2005 Jan 2;57(1):153-67. Review. PMID: 15518927
36. Impact of the new Beckman Coulter Cytomics FC 500 5-color flow cytometer on a regional flow cytometry clinical laboratory service. Luider J, Cyfra M, Johnson P, Auer I. *Lab Hematol.* 2004;10(2):102-8. PMID: 15224766
37. A combined in vitro/bioinformatic investigation of redox regulatory mechanisms governing cell cycle progression. Conour JE, Graham WV, Gaskins HR. *Physiol Genomics.* 2004 Jul 8;18(2):196-205. PMID: 15138307
38. Microtubule stress modifies intra-nuclear location of Msh2 in mouse embryonic fibroblasts. Marquez N, Chappell SC, Sansom OJ, Clarke AR, Teesdale-Spittle P, Errington RJ, Smith PJ. *Cell Cycle.* 2004 May;3(5):662-71. Epub 2004 May 15. PMID: 15044851
39. DRAQ5-based DNA content analysis of hematolymphoid cell subpopulations discriminated by surface antigens and light scatter properties. Yuan CM, Douglas-Nikitin VK, Ahrens KP, Luchetta GR, Braylan RC, Yang L. *Cytometry B Clin Cytom.* 2004 Mar;58(1):47-52. PMID: 14994375
40. Decreased uptake of bodipy-labelled compounds in the presence of the nuclear stain, DRAQ5. Snyder DS, Garon CF. *J Microsc.* 2003 Sep;211(Pt 3):208-11. PMID: 12950469
41. Optimization of three- and four-color multiparameter DNA analysis in lymphoma specimens. Plander M, Brockhoff G, Barlage S, Schwarz S, Rothe G, Kruechel R. *Cytometry A.* 2003 Jul;54(1):66-74. PMID: 12820122
42. Discrimination of DNA and RNA in cells by a vital fluorescent probe: lifetime imaging of SYTO13 in healthy and apoptotic cells. van Zandvoort MA, de Grauw CJ, Gerritsen HC, Broers JL, oude Egbrink MG, Ramaekers FC, Slaaf DW. *Cytometry.* 2002 Apr 1;47(4):226-35. PMID: 11933012
43. Characteristics of a novel deep red/infrared fluorescent cell-permeant DNA probe, DRAQ5, in intact human cells analyzed by flow cytometry, confocal and multiphoton microscopy. Smith PJ, Blunt N, Wiltshire M, Hoy T, Teesdale-Spittle P, Craven MR, Watson JV, Amos WB, Errington RJ, Patterson LH. *Cytometry.* 2000 Aug 1;40(4):280-91. PMID: 10918279
44. A novel cell permeant and far red-fluorescing DNA probe, DRAQ5, for blood cell discrimination by flow cytometry. Smith PJ, Wiltshire M, Davies S, Patterson LH, Hoy T. *J Immunol Methods.* 1999 Oct 29;229(1-2):131-9. PMID: 10556697