

UNIQUE Ageladine A

Novel non-toxic, cell permeable, pH-dependent Fluorescent Dye for Live Imaging

Unique Features of Ageladine A

- Wide range (pH 4 to pH 8). pH-dependent fluorescent dye in the blue-green range upon excitation with UV light. Stronger under acidic conditions and barely detectable in alkaline solutions.
- Non-toxic, highly cell/membrane permeable dye. No AM-esters or esterases involved. Trapped within the cells and acidic organelles through hydrophobic interactions with the inner side of the membranes.
- Barely metabolized and exerts long-term stability.
- Allows long term pH monitoring in acidic organelles, vesicles, cells, tissue and small animals over several days without side effects.

LITERATURE:

- Ageladine A, a pyrrole-imidazole alkaloid from marine sponges, is a pH sensitive membrane permeable dye: U. Bickmeyer, et al.; *BBRC* **373**, 419 (2008)
- Tracking of fast moving neuronal vesicles with ageladine A: U. Bickmeyer, et al.; *BBRC* **402**, 489 (2010)
- The alkaloid Ageladine A, originally isolated from marine sponges, used for pH-sensitive imaging of transparent marine animals: U. Bickmeyer; *Mar. Drugs* **10**, 223 (2012)
- Incorporated nematocysts in *Aeolidiella stephanieae* (Gastropoda, Opisthobranchia, Aeolidioidea) mature by acidification shown by the pH sensitive fluorescing alkaloid Ageladine A: D. Obermann, et al.; *Toxicon* **60**, 1108 (2012)

Ageladine A (synthetic)

AG-CMA-1001-C200

200 µg

Formula: C₁₂H₈Br₂F₃N₅O₂ · C₂HF₃O₂


MW: 357.0 · 114.0

CAS: 643020-13-7

Purity: ≥98%

Ex: between 325 & 415nm;
max. at 370nm

Em: from 415nm to >500nm;
max. at 415nm



APPLICATIONS:

Fluorescence microscopic monitoring of acidic organelles like lysosomes and endosomes, whole animals like plathelminths, cnidaria, larvae and eggs from different species • Screenings • Viability tests • Flow cytometry • Fluorescence lifetime imaging microscopy (FLIM)

Licensed from Alfred Wegener Institute (EP2156193A1; US8198098B2).




FIGURE:

CFLSM scan of *M. lignano* stained with Ageladine A (fluorescence is shown in false color). Ageladine A shows a pH-depending fluorescence. Highly fluorescent areas (white) are more acidic than dark areas (dark red).

See Backcover for Biological Activities of Ageladine A

Imaging Specifications

Loading times:

Cells in culture: ~10-30 minutes. Small animals (e.g. larvae, plathelminthes): ~30-120 minutes. Serum, salt content and ion assembly of the culture medium plays no role. All common buffers and culture media buffered up to pH 7.5 can be used.

Spectral properties:

Excitation between 325 & 415nm; max. at 370nm
Emission from 415nm to >500nm; max. at 415nm.

Experiments can be done with common filter settings (like for FURA-2). The fluorescence rises with lower pH in a linear range between pH 4 and pH 8. Quantitative results are achieved with ratiometric methods or FLIM.

Photobleaching occurs, but plays no role for monitoring at low excitation intensities and exposition times under usual laboratory conditions (even for several days).


Leakage was not observed after several days of incubation.

Toxic effects were not observed at concentrations up to 30µM. Patch clamp experiments with PC12-cells showed weak changes of the membrane potential at concentrations >10µM.


Recommended concentrations:

Between 1µM (cells) and 30µM (whole animals).

Imaging Examples



PC12-cells staining with ageladine A during UV excitation and the transmission image of the cells.



Nematocysts in different regions of *Aiptasia sp.* stained with Ageladine A: left: in acontia; right: in tentacles. Monitored by a confocal laser scanning microscope. Fluorescence is displayed in red (false colour).

Antiangiogenic Agent

BULK
available

Ageladine A is a unique brominated pyrrole-imidazole alkaloid, first isolated from the marine sponge *Agelas nakamurai*. It exhibits *in vitro* and *in vivo* **antiangiogenic activity**, which was initially considered to be associated with its moderate inhibition of various subtypes of matrix metalloproteinases (MMPs), but subsequently confirmed to result from its selective inhibition of kinases including dual specificity tyrosine-phosphorylation-regulated kinase (DYRK)1A, DYRK2, tyrosine kinase 2 (TYK2) and yeast Sps1/Ste20-related kinase 4 (YSK4). It is a highly selective **angiogenesis inhibitor** with no cytotoxicity against a panel of human cancer cell lines.



LITERATURE:

- Ageladine A: An antiangiogenic matrix metalloproteinase inhibitor from the marine sponge *Agelas nakamurai*: M. Fujita, et al.; *JACS* **125**, 1570 (2003)
- A one-pot synthesis and biological activity of ageladine A and analogues: S.R. Shengule, et al.; *Med. Chem.* **54**, 2492 (2011)
- Marine-derived angiogenesis inhibitors for cancer therapy: Y.Q. Wang & Z.H. Miao; *Mar. Drugs* **11**, 903 (2013)