

HYDROGEN PEROXIDE DETECTION

Red Hydrogen Peroxide Assay Kit

ENZ-51004


500 Reactions

HIGHLIGHT

- **Flexible:** Quantify hydrogen peroxide in solution, in cell extracts or directly from certain cells
- **Versatile:** Detect a variety of oxidase activities through enzyme-coupled reactions
- **Sensitive:** Detect as little as 10 picomoles of hydrogen peroxide
- **Complete:** Mix and read reagents, with minimal hands-on time
- **Convenient:** Homogenous assay format, fully compatible with HTS automation

Hydrogen peroxide (H_2O_2) is a reactive oxygen metabolic by-product that serves as a key regulator for a number of oxidative stress-related states. It is involved in a number of biological events and intracellular pathways that have been linked to several diseases. Measurements of reactive oxygen species help to determine how oxidative stress modulates varied intracellular pathways. The Red Hydrogen Peroxide Assay Kit uses a non-fluorescent peroxidase substrate to quantify hydrogen peroxide in solution, in cell extracts and released from certain live cells, such as granulocytes (Figure 1). It can also be used to detect a variety of oxidase activities through enzyme-coupled reactions. The kit provides an optimized and homogenous assay format that is compatible with automated liquid handling systems. The Red Hydrogen Peroxide Assay Kit offers a sensitive, one-step fluorometric assay to detect as little as 10 picomoles of H_2O_2 in a 100 μ L assay volume (Figure 2). The assay can be performed in a convenient 96-well or 384-well microplate format and is amenable to either fluorescence or absorbance based detection platforms with Ex/Em = 570 nm/590 nm or absorbance at 570 nm.

A.



B.




FIGURE 1: Conversion of Red Peroxidase Substrate into Resorufin: Horseradish peroxidase uses Red Peroxidase substrate as an electron donor during the reduction of hydrogen peroxide to water. The resultant product, resorufin, is a highly colored and fluorescent compound. Use of assay for detection of glucose oxidase activity (A). Use of assay in ELISA to detect HRP-antibody conjugates (B).




FIGURE 2: H₂O₂ dose response in a 384-well black well plate using a fluorescence microplate reader (NovoStar Labtech) measured with the Red Hydrogen Peroxide Assay Kit. As little as 0.01 µM H₂O₂ can be detected within a 30 minute incubation time (n=3). The inset shows low levels of H₂O₂ detection (< 1 µM).

Related Products

Product	Prod. No.	Size
ROS/RNS Detection Kit	ENZ-51001	200 Reactions
Luminol [3-Aminophthalhydrazide]	ENZ-52354	1 g

Switzerland & Rest of Europe

ENZO LIFE SCIENCES AG

Industriestrasse 17, Postfach
CH-4415 Lausen / Switzerland
Tel. + 41/0 61 926 89 89
Fax + 41/0 61 926 89 79
info-ch@enzolifesciences.com

North/South America

ENZO LIFE SCIENCES INTERNATIONAL, INC.

5120 Butler Pike
Plymouth Meeting, PA 19462-1202
USA
Tel. 1-800-942-0430 / (610) 941-0430
Fax (610) 941-9252
info-usa@enzolifesciences.com

Benelux

ENZO LIFE SCIENCES BVBA

Melkerijweg 3
BE-2240 Zandhoven / Belgium
Tel. +32/0 3 466 04 20
Fax +32/0 3 466 04 29
info-be@enzolifesciences.com

Germany

ENZO LIFE SCIENCES GmbH

Marie-Curie-Strasse 8
DE-79539 Lörrach / Germany
Tel. +49/0 7621 5500 526
Toll Free: 0800 6649518
Fax +49/0 7621 5500 527
info-de@enzolifesciences.com

UK & Ireland

ENZO LIFE SCIENCES (UK) LTD.

Palatine House
Matford Court
Exeter EX2 8NL / UK
Tel. 0845 601 1488 (UK customers)
Tel. +44/0 1392 825900 (overseas)
Fax +44/0 1392 825910
info-uk@enzolifesciences.com

For Local Distributors please visit our Website.