

# CUT&Tag-IT™ Assay Kit

Rapid and robust genome-wide analysis of histone marks at lower sequencing depths



**Cleavage Under Targets and Tagmentation (CUT&Tag)** is a method to investigate genomic localization of histone modifications and some transcription factors that reveals interactions between proteins and DNA or identifies DNA binding sites for proteins of interest.

Unlike MNase-Seq or ATAC-Seq methods that target open chromatin and are therefore dependent on chromatin accessibility, CUT&Tag utilizes an antibody-based enzyme tethering strategy to target specific histone modifications or proteins to reveal chromatin-binding information that is specific to those sites or proteins of interest.

CUT&Tag is based on the same principles as ChIP-Seq, but with several changes to the protocol that are advantageous in certain situations. Instead of the sonication of fixed chromatin and immunoprecipitation steps performed in ChIP-Seq protocols, in CUT&Tag, fresh (not frozen) unfixed cells are bound to concanavalin A beads and the antibody incubation is performed with cells in their native state. Directly following antibody binding, the chromatin is digested and NGS libraries are prepared in a single step by tagmentation using the protein A-Tn5 (pA-Tn5) transposome enzyme that has been pre-loaded with sequencing adapters.

CUT&Tag can rapidly produce high-quality results from less starting material than ChIP-Seq, and enables robust analysis from lower sequencing depths, saving both time and money.

## CUT&Tag vs. CUT&RUN vs. ChIP-Seq

	CUT&Tag-IT™ Assay Kit	CUT&RUN	ChIP-Seq
<b>Performed Under Native Conditions?</b>	Yes	Yes	No
<b>Chromatin Fragmentation Method</b>	Tn5-based tagmentation	MNase digestion	Sonication
<b>Cell Number Requirements</b>	5,000-500,000 cells	500,000 cells	1-10 million cells
<b>Sequencing Depth Required *</b>	2 million reads	8 million reads	20-50 million reads
<b>Integrated Library Preparation?</b>	Yes, uses tagmentation	No, separate library prep required	No, separate library prep required
<b>Compatible Targets</b>	Primarily histone modifications, some transcription factors and co-factors	Wide range of histone modifications, transcription factors, and co-factors	Wide range of histone modifications, transcription factors, and co-factors
<b>Workflow Length</b>	1-2 days	1-2 days	2-3 days

\* Kaya-Okur et al. *Nature Communications* (2019) 10:1930

## CUT&Tag-IT™ Assay Kit Advantages:

- Compatible with as few as 5,000 cells
- Complete kit with optimized protocol
- Developed for histone marks and some transcription factors
- Sequencing-ready libraries without the laborious and costly steps of ChIP-Seq
- Low background signal enables lower sequencing depth
- No artifacts caused by formaldehyde crosslinking



CUT&Tag-IT™ Assay Kit, Anti-Rabbit 16 rxns 53160


CUT&Tag-IT™ Assay Kit, Anti-Mouse 16 rxns 53165

### Don't forget to use an antibody validated for CUT&Tag!

Active Motif specializes in manufacturing high-quality antibodies to histones, histone modifications, chromatin proteins, and other factors, including a growing list of antibodies that we have experimentally validated in-house to work well in CUT&Tag assays. [Check out our GROWING list of CUT&Tag-validated antibodies:](#)

- AbFlex® CTCF antibody (rAb)
- CTCF antibody (pAb)
- Histone H3K4me1 antibody (pAb)
- Histone H3K4me2 antibody (pAb)
- Histone H3K4me3 antibody (pAb)
- Histone H3K9ac antibody (pAb)
- Histone H3K9me3 antibody (pAb)
- Histone H3K9me3 antibody (pAb)
- AbFlex® Histone H3K27ac antibody (rAb)
- Histone H3K27ac antibody (pAb)
- Histone H3K27ac antibody (pAb)
- Histone H3K27me3 antibody (pAb)
- Histone H3K36me3 antibody (pAb)
- AbFlex® RNA Pol II antibody (rAb)

## How the CUT&Tag-IT™ Assay Works



[www.activemotif.com](http://www.activemotif.com)



**Vinci-Biochem Srl**  
Via Ponte di Bagnolo, 10  
50059 Vinci (Firenze) Italy  
Tel: 0571 568 147 e 0571 568 135  
[vb@vincibiochem.it](mailto:vb@vincibiochem.it) – <https://www.vincibiochem.it>